

Biochemistry Analyzer



Product Specification

Automation	Semi Automatic
Brand	Sinothinker
Parameters	Diabetes Profile, Kidney Profile, Liver Profile, Lipid Profile, Cardiac Profile and Specific Protein
Assays	Clinical Chemistry
Power Supply	150V AC
Reaction Cuvette	Flow Cell, Reusable Cuvette
Test Method	End Point, Kinetic, Fixed Time and Absorbance
Minimum Order Quantity	2



ACCUCARE
We Care for Lives

Product Description

Product Details:

7 inch Color LCD, Touch screen

10 incubators(37±0.3 3 timers)

Test method Kinetic, Endpoint, Twopoint,Absorbance

Abs range:0~4.500Abs;Real time curve showing

Flow cell and Cuvette test mode compatible

7 wavelengths filter: 340,405,492,510,546,578,630nm

80 test item pre-programmed

Memory for 10,000 sample results

Light source hibernate function making lamp life longer

RS232 interface, PC connecting



We Care for Lives

ACCUCARE DIAGNOSTICS

Address: PLOT NO. 86, near PURI DIPLOMATIC GREENS,
Sai Kunj, Block A, New Palam Vihar Phase 3, Sector 110A,
Gurugram, Haryana 122017

www.accucare.in

